








1. Open the file of IVS-KA6000 EXE Install Main Program Installation.



2. Select "Yes(是)" to agree to the system installation program.




3. The instructions will guide the installation, click "Next" until the installation program starts.







4. Click "Finish" button to complete the installation process.



5. The shortcut of the IVS-KA6000 software will be displayed on the desktop, double click to run it.



6. Since the IVS-KA6000 software will generate parameters, your operating system may pop up a safety warning. IVS-KA6000 will not have any impact on the computer system, please select "Yes(是)".





7. Start to use IVS-KA6000. Insert the USB KEY into the computer to use it, or press "Free Trial" to try it for 14 days

